



technical documentation

GANSER
LIFTSYSTEME

GTL 30



Ganser Maschinen GmbH
Markt 26
4171 St. Peter a. Wbg.
AUSTRIA



office@ganserlifte.at



0043 7282 80 71



Table of contents

1.	Ganser lifting systems – From vision to reality.....	2
2.	Advantages of the GTL30.....	2
3.	Stairlift design and operating control options	3
3.1.	Operating control options	3
3.2.	Parts of platform stairlift.....	4
4.	Accessories.....	5
4.1.	Emergency telephone	5
4.2.	Permanently mounted external control panels.....	5
4.3.	Handheld remote control	5
4.4.	Cable remote.....	5
4.5.	Front access ramp	6
4.6.	Folding seat	6
5.	External controls /variants	7
5.1.	Surface mounted (Standard)	7
5.2.	Flush mounted	7
6.	General technical data	8
7.	Installations	9
7.1	Wall mounted / Examples	9
7.2	Support assembly / Examples	10
8	Sample drawings.....	11
8.1	Wall mounted / Sample drawing	11
8.2	Support assembly / Sample drawing	12
9	Minimum dimensions for various platform sizes	13
9.1	180° curve with inclination	13
9.2	180° curve and wall bevelled (in the area from the rail)	14
9.3	180° curve with longer lugs	15
9.4	180° curve with columns	16
9.5	90° curve with longer lugs	17
9.6	90° curve and wall bevelled (in the area from the rail)	18
9.7	90° curve and columns.....	19
10	Stress points.....	20
11	Electrical drawing.....	21
12	Type examination	25

1. Ganser lifting systems – From vision to reality

The success of Ganser Liftsysteme as a stairlift manufacturer is rooted in many years of company development that goes back to 1965. Over the years, the company has positioned itself on the market as the most important partner in the area of lift development and lift production. This goes hand in hand with the vision of enabling people all over the world to have a barrier-free life. Together with customer-oriented company management, excellent employees and a global dealer network, we have already made the vision of barrier-free living in many homes, public institutions and companies a reality. The team at Ganser Liftsysteme stands for flexibility, expertise and reliability and thus enables customers worldwide to benefit from individually adapted lift systems and to live barrier-free. Ganser Liftsysteme also turns your vision into reality - entrust us with your challenges.

2. Advantages of the GTL30

- Upper stainless steel tube can be used as a handrail
- Minimum installation depth of 350 mm in the park position
- EN 81-40
- The most difficult roadways possible, including negative curves
- High quality weatherproof materials
- New control including display for quick error detection
- Soft start at the beginning, speed reduction in the curves and when changing the incline
- Suitable for narrow stairs

3. Stairlift design and operating control options

3.1. Operating control options



Figure 1: Control panel on lift



Figure 2: Permanently mounted external remote control panel



Figure 3: Handheld remote control



Figure 4: Additional cable remote for folding seat

01	DOWN control key	Depending on the direction of the arrow, the stairlift moves up or down.
02	UP control key	
03	Open platform	Depending on the direction of the arrow, the platform will be folded up or down.
04	Close platform	
05	Emergency alarm	Activates an emergency alarm sound
06	Emergency stop button	When activated, the stair lift stops. The button must be turned to release.
07	Key switch	The system can only be operated with the key.

Table 1: List of operating controls

3.2. Parts of platform stairlift



Figure 5: Parts of stairlift

08	Handrail
09	Safety barrier
10	Access ramp
11	Roll-off barrier, front access ramp (optional)
12	Bottom contact sensor
13	Intelligent display
14	Control panel on lift
15	Service access (rear of lift equipment)
16	Side safety cut-out switch
17	Switch to power on stairlift.
18	Locking safety gear (rear of lift equipment)

Table 2: Parts of platform stair lift

4. Accessories

Ganser stairlifts can also be equipped with optional accessories to better meet customer requirements.

4.1. Emergency telephone

The stairlift can be equipped with an optional GSM telephone, so that there is always a telephone on the lift. We recommend that every user of this stairlift system, especially in the private sector, carries a mobile telephone whilst operating the system.



4.2. Permanently mounted external control panels

External control panels (see Figure 2) can be activated or deactivated with a key. Use of this key ensures that unauthorised persons cannot use the stairlift. External control panels are fixed to the wall at each station.



4.3. Handheld remote control

This device allows the operator as much free space as possible while operating the stairlift. They will no longer be dependent on the stationary control panels at the stations or on the carriage but can control the lift with a small remote control (see Figure 3) which they can always carry on their person. The receiver is mounted on the carriage and transmits the radio commands to the control unit.

Due to the absence of external control panels, the system is largely protected from vandalism and deliberate damage.



4.4. Cable remote

The platform stair lift can optionally be equipped with an additional cable remote control (see Figure 4), which ensures ergonomic operation when using the folding seat.



4.5. Front access ramp

Sometimes space limitations in front of the first step mean it is not possible to drive onto the platform using the two standard access ramps. In these cases, an additional access ramp must be fitted along the long edge of the platform, with a roll bar on the lower side of the platform if required.




4.6. Folding seat

The folding seat allows the operator to sit while travelling on the stairlift. It is designed with a comfortable plastic surface. When the seat is not in use, it can be folded up to save space.




5. External controls /variants

5.1. Surface mounted (Standard)



5.2. Flush mounted



6. General technical data

Electrical supply	230V / 50 Hz
Number of trips	approx. 10 / hour
Motor power	0.6 kW
Nominal current	max. 13 A (depending on the version)
Starting current	16 A (depending on the version)
Speed	approx. 0.15 m / sec (soft start at start / stop)
Drive	rack with gear
Variant	battery version
Roadway	stainless steel tube 42.2 mm, can also be used as a handrail
Slope	up to 60°
Load capacity	up to 300 kg
Platform size	Standard platform size 1000 x 800 mm (smaller sizes at no extra charge), max. 1000 x 1000 mm (special sizes also possible)
Color	all RAL colors possible, special colors on request or stainless steel in K400 or K800

7. Installations

7.1 Wall mounted / Examples




7.2 Support assembly / Examples




8 Sample drawings

8.1 Wall mounted / Sample drawing




8.2 Support assembly / Sample drawing



9 Minimum dimensions for various platform sizes

9.1 180° curve with inclination


	600x900	650x900	700x900	750x900	800x900	800x1000	800x1200
A: platform closed [mm]	395	395	395	395	395	395	395
B: Plattform opend [mm] *	886	936	986	1036	1086	1086	1086
C: minimum staircase [mm] *	1020	1060	1110	1150	1195	1220	1260
D: minimum slewing range [mm] *	1160	1200	1250	1290	1335	1360	1400



* with front ramp +100mm

9.2 180° curve and wall bevelled (in the area from the rail)


	600x900	650x900	700x900	750x900	800x900	800x1000	800x1200
A: platform closed [mm]	395	395	395	395	395	395	395
B: Plattform opend [mm] *	886	936	986	1036	1086	1086	1086
C: minimum staircase [mm] *	1020	1060	1110	1150	1195	1220	1260
D: minimum slewing range [mm] *	1020	1060	1110	1150	1195	1220	1260



* with front ramp +100mm

9.3 180° curve with longer lugs


	600x900	650x900	700x900	750x900	800x900	800x1000	800x1200
A: platform closed [mm]	435	435	435	435	435	435	435
B: Plattform opened [mm] *	1026	1076	1126	1176	1226	1226	1226
C: minimum staircase [mm] *	1060	1100	1150	1190	1235	1260	1300
D: minimum slewing range [mm] *	1080	1120	1170	1210	1255	1280	1320



* with front ramp +100mm

9.4 180° curve with columns


	600x900	650x900	700x900	750x900	800x900	800x1000	800x1200
A: platform closed [mm]	449	449	449	449	449	449	449
B: Plattform opened [mm] *	940	990	1040	1090	1140	1140	1140
C: minimum staircase [mm] *	1080	1120	1170	1210	1255	1280	1320
D: minimum slewing range [mm] *	1080	1120	1170	1210	1255	1280	1320



* with front ramp +100mm

9.5 90° curve with longer lugs


	600x900	650x900	700x900	750x900	800x900	800x1000	800x1200
A: platform closed [mm]	415	415	415	415	415	415	415
B: Plattform opened [mm] *	906	956	1006	1056	1106	1106	1106
C: minimum staircase [mm] *	1040	1080	1130	1170	1215	1240	1280
D: minimum slewing range [mm] *	1040	1080	1130	1170	1215	1240	1280



* with front ramp +100mm

9.6 90° curve and wall bevelled (in the area from the rail)


	600x900	650x900	700x900	750x900	800x900	800x1000	800x1200
A: platform closed [mm]	395	395	395	395	395	395	395
B: Plattform opend [mm] *	886	936	986	1036	1086	1086	1086
C: minimum staircase [mm] *	1020	1060	1110	1150	1195	1220	1260
D: minimum slewing range [mm] *	1020	1060	1110	1150	1195	1220	1260



* with front ramp +100mm

9.7 90° curve and columns

	600x900	650x900	700x900	750x900	800x900	800x1000	800x1200
A: platform closed [mm]	449	449	449	449	449	449	449
B: Plattform opened [mm] *	940	990	1040	1090	1140	1140	1140
C: minimum staircase [mm] *	1080	1120	1170	1210	1255	1280	1320
D: minimum slewing range [mm] *	1080	1120	1170	1210	1255	1280	1320



* with front ramp +100mm

10 Stress points

forces GTL30

Total weight [kg]	170
load capacity [kg]	200
a [mm]	520
Fges [statisch] [N]	4630
Fh [N]	4363
Fv [N]	4630
Fs [N]	2160A

ATTENTION!
If this distance becomes larger, the forces are also increasing!

Trest Mounted


Wall Mounted

Mounting Methods:


- The scissor lift can be mounted on support pillars or directly on the wall.
- A mounting method should be chosen such that the applied loads are adequately supported.
- Support pillars can be arranged in a variety of ways to suit the application.
 - With footplate fixed directly to trestle.
 - With footplate fixed to trestle and additionally braced horizontally into wall.

		Auftrags-Nr.: 1881121101 01/18	
GANSENER LIFTSYSTEME WERKSTOFFE GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME		Material: 1 : 10	Masse:
GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME		Werkstoff:	Stück: 1
GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME		Nachbearbeitung:	GTL30
GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME		Refertigung GTL30	
GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME GANSENER LIFTSYSTEME		1	


11 Electrical drawing




control station on board




GANSEER MASCHINEN GmbH Central Office: Badenauweg 79 5-4131 Otterfingen/Donau Tel.: +49 7362 8071-0 info@ganseerlift.de www.ganseerlift.de		
G1130-V2-6N Serial no. 3-2023 page 2	G1130-V2-6N Serial no. 3-2023 page 2	G1130-V2-6N Serial no. 3-2023 page 2



external control station







GANSEK Maschinen GmbH
 Tempelwiese 11 · Heimbach
 A-4171 St. Peter/Obg. Markt 26
 Tel.: (+3 752 8771-0
 info@gansek.at · www.gansek.at

DT130-V2-eN
 main drive

5.7023	Seite
	3



charging



- 81 main fuse 60A
- 82 controller fuse 1A
- 83 motor fuse 5A
- 84 charging fuse resettable 2A
- 85 buzzer
- 86 flashing light (warning)
- 87 main contactor
- 88 safety contactor
- 89 main contactor
- 90 safety contactor
- 91 photo barrier
- 92 photo barrier
- 93 photo barrier
- 94 photo barrier
- 95 photo barrier
- 96 photo barrier
- 97 photo barrier
- 98 photo barrier
- 99 photo barrier
- 100 photo barrier
- 101 photo barrier
- 102 photo barrier
- 103 photo barrier
- 104 photo barrier
- 105 photo barrier
- 106 photo barrier
- 107 photo barrier
- 108 photo barrier
- 109 photo barrier
- 110 photo barrier
- 111 photo barrier
- 112 photo barrier
- 113 photo barrier
- 114 photo barrier
- 115 photo barrier
- 116 photo barrier
- 117 photo barrier
- 118 photo barrier
- 119 photo barrier
- 120 photo barrier
- 121 photo barrier
- 122 photo barrier
- 123 photo barrier
- 124 photo barrier
- 125 photo barrier
- 126 photo barrier
- 127 photo barrier
- 128 photo barrier
- 129 photo barrier
- 130 photo barrier
- 131 photo barrier
- 132 photo barrier
- 133 photo barrier
- 134 photo barrier
- 135 photo barrier
- 136 photo barrier
- 137 photo barrier
- 138 photo barrier
- 139 photo barrier
- 140 photo barrier
- 141 photo barrier
- 142 photo barrier
- 143 photo barrier
- 144 photo barrier
- 145 photo barrier
- 146 photo barrier
- 147 photo barrier
- 148 photo barrier
- 149 photo barrier
- 150 photo barrier
- 151 photo barrier
- 152 photo barrier
- 153 photo barrier
- 154 photo barrier
- 155 photo barrier
- 156 photo barrier
- 157 photo barrier
- 158 photo barrier
- 159 photo barrier
- 160 photo barrier
- 161 photo barrier
- 162 photo barrier
- 163 photo barrier
- 164 photo barrier
- 165 photo barrier
- 166 photo barrier
- 167 photo barrier
- 168 photo barrier
- 169 photo barrier
- 170 photo barrier
- 171 photo barrier
- 172 photo barrier
- 173 photo barrier
- 174 photo barrier
- 175 photo barrier
- 176 photo barrier
- 177 photo barrier
- 178 photo barrier
- 179 photo barrier
- 180 photo barrier
- 181 photo barrier
- 182 photo barrier
- 183 photo barrier
- 184 photo barrier
- 185 photo barrier
- 186 photo barrier
- 187 photo barrier
- 188 photo barrier
- 189 photo barrier
- 190 photo barrier
- 191 photo barrier
- 192 photo barrier
- 193 photo barrier
- 194 photo barrier
- 195 photo barrier
- 196 photo barrier
- 197 photo barrier
- 198 photo barrier
- 199 photo barrier
- 200 photo barrier
- 201 photo barrier
- 202 photo barrier

Ganser Maschinen GmbH
Zweigstelle u. Labordosen
A-4150 St. Pölten/Donau Markt 26
Tel.: +43 7282 8071-0
Telefax: +43 7282 8071-20
E-Mail: info@ganser.com

Legend	4.0028
GT130-V2-80	Schüler
Seite 4	1
1	2

12 Type examination



EG-Baumusterprüfbescheinigung Certificate of EC type examination

Nr.: TÜV-A-MHF/MG 09-03488 VV

gemäß Richtlinie 2006/42/EG (Maschinen) - Verlängerung
according to the Directive 2006/42/EC (Machinery) - extension

Auftraggeber: Applicant:	Ganser Maschinen GmbH Markt 26 A-4171 St. Peter am Wimberg	Hersteller: Manufacturer:	siehe Auftraggeber see applicant
Produkt: Product:	Treppenschrägaufzug für den Personentransport	Zubehör: Accessories:	--
Typ: Type:	Plattform-Treppenlift GTL-30, GTL 30 Akku		
Beschreibung: Description:	Plattform-Treppenlift mit Kurvenfahrt für Personen mit eingeschränkter Mobilität, mit integrierter Fangvorrichtung, Geschwindigkeit < 0,15 m/sek., Tipfbetrieb		

Inspektionsgrundlage gemäß EN ISO/IEC 17020: Richtlinie 2006/42/EG idgF
Tested according to Inspection Body EN ISO/IEC 17020: Directive 2006/42/EC in the current version

Mitgeltende Prüfgrundlagen: EN ISO 12100, EN ISO 13849, EN ISO 13850, EN ISO 13857,
EN ISO 14118, EN ISO 13854, EN 60204-1, EN 81-40
Applicable standards: (Anwendung gemäß Risikobeurteilung / application according to risk assessment)

Bemerkungen: Hebeanlage und Fangvorrichtung sind als geprüfte Einheit zu sehen
Remarks:

Hiermit bestätigt die TÜV AUSTRIA SERVICES GMBH als Notifizierte Stelle (ID-Nr. 0408), dass das oben angeführte Produkt den grundlegenden Sicherheits-Anforderungen der Richtlinie 2006/42/EG entspricht. Grundlage dieser Bescheinigung ist das zur Prüfung vorgelegte Prüfmuster und die technische Dokumentation.

Gemäß Artikel 5 und Anhang III ist am Produkt die CE-Kennzeichnung vorgesehen.

Hereby TÜV AUSTRIA SERVICES GMBH certifies as Notified Body (ID-No 0408), that the above-mentioned product meets the essential safety requirements of the Directive 2006/42/EC. The certificate is based on the test specimen and the technical documentation subjected to the test.

According to Article 5 and annex III the CE mark is foreseen to be affixed on the product.

Prüfbericht: TÜV-A/MHF/MG 09-03488 v. 05.01.2010
Test Report:

26.01.2010 / 02.04.2015
Datum Ausstellung / Verlängerung V
Date of Issue / extension V

01.02.2020
Datum
date



31.01.2025
verlängert bis
extended till

Auszugsweise Vervielfältigung nur mit Genehmigung der TÜV AUSTRIA SERVICES GMBH gestattet
Duplication of this document in parts is subject to the approval by TÜV AUSTRIA SERVICES GMBH

09-03488VV cba GANSENER GTL-30-MRL-MG-SMF_Verl.docx
FM-NE-AS-MRL-0100_Verlängerung
Rev 03
Seite 1 von 3

TÜV AUSTRIA SERVICES GMBH
Industry and Energy
Anlagen- und Maschinensicherheit

Deutschstraße 12
A-1230 Wien
Tel.: +43 (0) 464 0203
http://www.tuv.at

ZERTIFIKAT | CERTIFICATE | CERTIFICAT | CERTIFICADO | СЕРТИФИКАТ | 證書 | 인증서

Gültigkeit von EG-Baumusterprüfbescheinigungen

1. Erlöschen von Bescheinigungen

Eine Bescheinigung wird auch ohne besondere Mitteilung der TÜV AUSTRIA SERVICES GMBH ungültig, wenn die auf der Bescheinigung angegebene Gültigkeitsdauer abgelaufen ist.

2. Einschränkung, Aussetzung, Ungültigerklärung, Rückzug von Bescheinigungen

Die Bescheinigungen können von der Notifizierten Stelle (TÜV AUSTRIA SERVICES GMBH) mit sofortiger Wirkung eingeschränkt, ausgesetzt oder für ungültig erklärt und zurückgezogen werden, wenn:

- (a) das bescheinigte Produkt nicht mehr dem genehmigten Muster entspricht,
- (b) Produkte für Endbenutzer oder Dritte eine Gefährdung darstellen,
- (c) zum Zeitpunkt der Prüfung Tatsachen nicht oder nicht richtig gesehen und beurteilt worden sind oder auch nicht erkennbar waren, die einer positiven Bewertung entgegengestanden hätten - hierzu gehört z. B. auch eine fehlerhafte Kategorisierung von Produkten in bestimmte Risikoklassen oder die Einordnung nach Verwendungszweckarten und zwar auch aufgrund eines von der Notifizierten Stelle (TÜV AUSTRIA SERVICES GMBH) zu verantwortenden Fehlers oder Mangels bei der Prüfung,
- (d) bei wiederkehrenden Überwachungen, bei Marktkontrollen oder sonst sich nachträglich herausstellenden Produkt- oder Systemmängel, die nicht vom Bescheinigungsinhaber in einer angemessenen Frist abgestellt werden,
- (e) der Bescheinigungsinhaber die wiederkehrenden Überwachungstätigkeiten der Notifizierten Stelle (TÜV AUSTRIA SERVICES GMBH) nicht durchführen lässt oder die ordnungsgemäße Durchführung behindert oder einschränkt,
- (f) Bescheinigungen oder Kopien von Bescheinigungen geändert und damit gefälscht worden sind,
- (g) irreführende oder anderweitig unzulässige Werbung mit den Bescheinigungen betrieben wird,
- (h) fällige Entgelte für die Bescheinigungen und/oder Produkt-Prüfung vom Bescheinigungsinhaber nicht in der gestellten Frist entrichtet werden.

3. Verfahren zum Bescheinigungs-Entzug

Die Bescheinigung ist im Falle der Ungültigkeitserklärung unverzüglich im Original an die TÜV AUSTRIA SERVICES GMBH zurückzugeben.

Der Notifizierten Stelle (TÜV AUSTRIA SERVICES GMBH) wird die Veröffentlichung von Einschränkungen, Aussetzungen, Ungültigkeitserklärungen und Zurückziehungen sowie Löschungen von Bescheinigungen vom Auftraggeber gestattet.

Validity of EC Type Examination Certificates

1. Expiry of certificates

A certificate expires even without special notice by TÜV AUSTRIA SERVICES GMBH when the expiry date specified on the product certificate is reached.

2. Restriction, suspension, cancellation, withdrawal of certificates

Certificates can be restricted, suspended or cancelled and withdrawn with immediate effect by the Notified Body (TÜV AUSTRIA SERVICES GMBH) if:

- (a) the certified product no longer corresponds to the approved type,
- (b) products represent a risk for end users or third parties,
- (c) facts precluding positive certification were not correctly seen or assessed at the time of the test or were not recognizable - this also includes e.g. incorrect categorization of products in defined risk classes or the classification according to types of intended uses, also if this is due to an error or fault of the certification for which the Notified Body (TÜV AUSTRIA SERVICES GMBH) is responsible,
- (d) defects of the product or of the system are detected in the course of repeated tests, market inspections or in any other way at a later moment and are not eliminated within a reasonable period by the holder of the certificate,
- (e) the holder of the certificate omits to have the repeated inspection carried out by the Notified Body (TÜV AUSTRIA SERVICES GMBH) or impedes or restricts its execution in due form,
- (f) certificates or copies of certificates were modified and therefore forged,
- (g) misleading or otherwise inadmissible promotion is made with product certificates,
- (h) due fees for the certification and/or product test are not paid by the holder of the certificate within the defined period.

3. Certificate withdrawal procedure

If the certificate is declared invalid, its original copy has to be returned to the Notified Body (TÜV AUSTRIA SERVICES GMBH).

The customer authorizes the Notified Body (TÜV AUSTRIA SERVICES GMBH) to publish restrictions, suspensions, cancellations and withdrawals as well as annulments of certificates.

Anhang zur EG-Baumusterprüfbescheinigung
Attachment to Certificate of EC-type examination
Nr.: TÜV-A-MHF/MG 09-03488 VV

Schematic representation of platform staircase lift GTL30

